

# Advanced in vitro exposure systems

VITROCELL® 12/6 CF Stainless Steel




# VITROCELL® 12/6 CF Stainless Steel Exposure Module

For 6 cell culture inserts (12-well size)

The VITROCELL® 12/6 CF module system has been specifically designed and engineered to facilitate the research of human cell cultures in direct exposure to airborne substances such as gases, complex mixtures, nanoparticles and fibers. The system authentically simulates the conditions of human physiological exposure.


Normally 3 compartments are used for exposure to the substances and 3 compartments for clean air control. The capacity can be easily increased by adding more modules. Each module can receive a different dose concentration, so that a complete dose/response profile is obtained in one experiment.

The cells are exposed at the air/liquid interface on 6 cell culture inserts using low flow rates of the aerosol. After exposure, the cells are further processed to measure a wide range of endpoints, e. g. cytotoxicity, genotoxicity, proliferation, cellular and oxidative stress as well as inflammation (see also VITROCELL® assay guide).



## Culture media supply


The media is supplied to the module by gravimetric method, syringe or using a media pump for intermittent/continuous media exchange.



VITROCELL® 12/6 CF Dynamic Dilution System and rack



Exposure Systems



Dose Supply



Dilution



Racks &amp; Carts



Dose Monitoring



Skin Exposure



Auxiliary Equipment

## Base module

The VITROCELL® 12/6 CF modules casing is made of electropolished stainless steel. It is designed with six compartments for 12-well cell culture inserts and is fully autoclavable at 121° C (250° F) for 20 min. Through the employment of a special adapter set, this module can also be flexibly combined for use with 24-well cell culture inserts. A constant unit temperature is guaranteed using a regulated flow of temperature-controlled water through the module. The media are supplied individually to each well compartment. Optionally, media exchange can be carried out on a continuous basis using a precision media pump.

This module is renowned in scientific circles for its reliability and durability.



Suitable for six cell culture insert (12-well size) – adaptor set for 24-well sized inserts available.

## Aerosol exposure top with individual aerosol inlets

Special sealing and connection mechanisms guarantee a hermetic connection of the base module with the aerosol exposure top.

The aerosol inlet-stream flows through specifically shaped inlets. They are available in stainless steel or stainless steel with unique VITROGLIDE surface treatment for lowest adhesion results. The stainless steel / VITROGLIDE inlets are specifically designed for aerosols containing nano-particles.

The aerosol inlets are connected to the distribution or dilution systems. Extraction from the module takes place via small holes using a controlled vacuum flow. VITROCELL® calibration valves, which are connected to a vacuum pump, control flow rates effectively. The aerosol exposure top features an in/outlet for heating the water circuit, ensuring corresponding temperatures with the base module. It is made of high-quality anodized aluminum.



The module can easily be locked hermetically with the two levers attached to the base.

## Option for microbalance sensor


The QCM sensor can be integrated in the VITROCELL® 12/6 CF exposure module. It is capable of measuring the deposited mass at a resolution of 10 nanogram/cm<sup>2</sup> per second. Results are reported online by the VITROCELL® Monitor software. Data is presented in graphs and stored in MS Excel®.




## System available in different setups for dosimetry

The system can be combined with a separate single module for dosimetry. The dosimetry module can accommodate a stainless steel insert for trapping constituents of the aerosol in liquids. Alternatively it can be equipped with a microbalance sensor.

A photometer can be added either directly in the aerosol flow of the Exposure Module or at a separate port of the main aerosol flow.



Rack with Dosimetry Module and Photometer connected to Exposure Module.



Rack with Dosimetry Module and Photometer connected to the main aerosol flow.


## Features

- Suitable for **COSTAR®**, **FALCON®** and **ThinCert®** 12-well sized cell culture inserts
- For higher throughput
- Autoclavability of all components
- Base module made of electropolished stainless steel for extreme durability
- Base module with water heating base plate
- Hyperboloid geometry of inlets for optimized particle deposition and distribution
- Option: Quartz Crystal Microbalance

# VITROCELL® Exposure Systems for Inhalation Toxicology


## Direct Exposure Technology at Air/Liquid Interface

In response to the scientific need to expose in physiologically relevant conditions, VITROCELL® exposure modules have been specifically designed and engineered to enable direct exposure of mammalian cells or tissue at the air/liquid interface where the cell systems are not covered with culture media. Researchers can thus use all cell types cultivated on microporous membranes. This approach allows for more credible and authentic results than by submerged exposure due to a closer replication of the human physiology.



### The advantages:

- No losses
- No dissolution
- No reaction of constituents with culture media
- High sensitivity



Optimized  
flow dynamics





## VITROCELL® EXPOSURE AT THE AIR/LIQUID INTERFACE

**Phases**

Gaseous  
Semi-volatile  
Particle

**A**  
**B**  
**C**

**Submerged Cultivation and Exposure in Incubator**

**A** Media above cells  
**B** Cells on membrane  
**C** Media below cells

Interaction of test components with culture media

**Low sensitivity**

**Phases**

Gaseous  
Semi-volatile  
Particle

**A**  
**B**  
**C**

**Suspension Cultivation and Exposure in Incubator**

**A** Cells in media  
**B** Cells on membrane  
**C** Media below cells

Interaction of test components with culture media

**Low sensitivity**

**Phases**

Gaseous  
Semi-volatile  
Particle

**A**  
**B**  
**C**

**Air / Liquid Cultivation and Exposure in Exposure Module**

**A** Direct and controlled exposure of test atmosphere to cells  
**B** Cells on membrane  
**C** Media below cells

No losses  
No reaction of principle components with culture media

**High sensitivity of system**

The exposure of mammalian cells or tissues to airborne substances is frequently performed under submerged conditions. In doing so, the test substances are dosed into the culture media. This procedure results in an undesired interaction of the formerly airborne substances with the media, causing limitations for authentic analysis.

Therefore VITROCELL® recommends the air/liquid interface exposure technology.




## About VITROCELL®

**VITROCELL® exclusively concentrates on the developing, producing, installing, training and servicing of advanced *in vitro* exposure systems.**

The VITROCELL® Systems' team is driven by their vision for new in-vitro standards through state-of-the-art technology, highly qualified workmanship and absolute client dedication. VITROCELL® has successfully collaborated with clients from leading research institutes, contract research organizations, regulatory authorities or industrial laboratories across the world. Working with our team experts, all modules have been tailored to create durable and complete turnkey-systems for *in vitro* inhalation toxicology. Gases, environmental atmospheres, nano particles and complex mixtures are analyzed on lung cells at the air/liquid interface using these systems. VITROCELL® technologies are also applicable to solutions for skin research.

Over a decade of devotion to research in this specific field has given our team of design & precision manufacturing specialists the opportunity to mentor highly diversified and complex projects **from conception to completion**. We strive to become a constructive member of each research team, providing support when it is needed, advice when it is required and modules of the highest quality, which are even polished by hand before leaving here to be integrated into your workspace. Every piece of our German engineered equipment is manufactured to the highest of standards – yours.

For more information  
please scan the QR-Code:



VITROCELL® Systems GmbH  
Fabrik Sonntag 3  
79183 Waldkirch  
Germany

Tel. +49 7681 497 79-50  
Fax +49 7681 497 79-79  
Email: [info@vitrocell.com](mailto:info@vitrocell.com)  
[www.vitrocell.com](http://www.vitrocell.com)

