

# VITROCELL® Cloud Alpha

A Family of Exposure Systems Suitable for Nebulized Solutions and Suspensions

Patrick Weindl, Heidi Ortolf-Wahl, Tobias Krebs, VITROCELL Systems GmbH, 79183 Waldkirch, Germany

## Optimal Solution for Everyday Experiments at the Air/Liquid Interface

The VITROCELL® Cloud Alpha series is our newest innovation and presents a great leap forward in automated exposure of cell cultures. It combines reliable exposure of cell cultures from the respiratory tract with ease of use. The series comprises the Cloud Alpha 6, -12, -96 and Cloud Alpha MAX – they represent an optimal solution for everyday experiments at the Air/Liquid Interface using 6-well, 12-well, 24-well or 96-sized inserts.

We have run extensive tests to ensure the Cloud Alpha series devices match our standards in terms of the established characteristics for Cloud principle exposure devices: Deposition Efficiency, Spatial Deposition and Repeatability.

## Methods

To determine the performance of the devices, Fluorescein-Sodium was nebulized with standard Aeronex vibrating mesh nebulizers, imitating the procedure of a regular cell culture experiment and following the well described Cloud principle of operation.<sup>[1]</sup> Previously reported optimal nebulization volumes (see tab. 1) were used. The Fluorescein was trapped in PBS using stainless steel inserts. The concentration of deposited fluorescein was then determined via spectrofluorimetry. This method has been described in detail before.<sup>[1][2]</sup>

The characteristic values may be calculated via:

$$(1) \text{ DepositionEfficiency} = \frac{C_{fluo.Insert} * V_{catch.Insert}}{C_{fluo.Nebulized} * V_{Nebulized} * \frac{A_{Insert}}{A_{Chamber}}}$$

$$(2) \text{ SpatialDeposition} = \sqrt{\frac{\sum_{i=1}^n (C_{fluo.Insert_i} - \mu_i)^2}{n}} * \frac{1}{\mu_i}$$

The shown efficiencies are values averaged over all inserts and over several experiments (minimum 3, maximum 9). The variance between this averaged efficiency per run has then been used as a measure for repeatability. To address the different ratios between exposed area and cell culture area per device, we propose the measure Insert Deposition Efficiency:

$$(3) \text{ InsertDepositionEfficiency} = \text{DepositionEfficiency} * \frac{A_{Culture}}{A_{Chamber}}$$



Fig. 1: VITROCELL® Cloud Alpha Family overview. Cloud Alpha 6, -12, -96+, -96 and Cloud Alpha MAX, left to right.




Fig. 2: Insert-to-insert variability per device type (n left to right: 4, 4, 4, 9, 7, 7, 6)




Fig. 3: Efficiency (overall and culture related) in available devices




Fig. 4: Experiment-to-experiment variability in available devices




Fig. 5: Cell culture efficiency

	Alpha 6		Alpha 12		Alpha 96+	Alpha 96	Alpha MAX	
	6-well	12-well	24-well	12-9				12-3
Deposition Efficiency	90%	84%	72%	77%	74%	76%	74%	32%
Exposure area / (cm <sup>2</sup> )	145	145	145	141,5	54	200,75	108,75	4,71
Culture growth area / (cm <sup>2</sup> )	28,0	6,7	2,0	10,1	3,4	13,7	13,7	4,7
Insert Deposition Efficiency	17%	4%	1%	5%	5%	5%	9%	32%
Optimal deposition for standard nebulization of 1% solution / (µg/cm <sup>2</sup> )	13,8	13,8	13,8	14,1	14,8	14,9	18,4	42,5
Nebulization volume / (µL)	min		50	50	30	300	30	10
	standard		200	200	80	300	200	20
	max		400	400	130	300	300	40

Tab. 1: Areas, Volumes and Efficiencies of the VITROCELL® Cloud Alpha devices

## Features:

- Cloud dynamics by substance-safe vibrating mesh nebulizers
- Programmable, automated exposures
- Very low insert-to-insert variability/high reproducibility
- Air Liquid Interface (ALI) cell culture exposure
- Direct guidance of test substance to the cells
- No external airflow required



Fig. 6: Touchscreen control of VITROCELL® Cloud Alpha devices




Fig. 7: Vibrating mesh nebulizers for aerosol generation

[1] Lenz et al, A dose-controlled system for air-liquid interface cell exposure and application to zinc oxide nanoparticles, Particle and Fibre Toxicology, page 3 ff., 2009 (DOI: 10.1186/1743-8977-6-32)

[2] Ding et al, Quartz crystal microbalances (QCM) are suitable for real-time dosimetry in nanotoxicological studies using VITROCELL® Cloud cell exposure systems, Particle and Fibre Toxicology, page 6 f., 2020 (DOI: 10.1186/12989-020-00376)