

Holder	Product Picture	Standard	Necessary Accessories	Products	Activation
Cigarettes		Yes	—	Conventional Cigarettes	Manually
Electronic cigarettes		No	Software Option (Square Profile)	Standard draw actuated E-Cigarettes	Draw actuated
ENDS cigarette (Electronic Nicotine Delivery)		No	Software Option (Square Profile) Seals depending on product	Differnt ENDS products possible (e.g. Juul, blu, Alto)	Draw actuated
Renova Zero (ENDS product)		No	Software Option (Square Profile)	Renova Zero product	Draw actuated
IQOS or GLO Heated Tobacco Product		No	Software Option (Square Profile)	IQOS or GLO Products	Manually
Cigar Holder System		No	—	Cigar like products OD 6-22 mm	Manually